

BRAF Mouse mAb

ABZ4240


Product Overview

Name	BRAF Mouse mAb
Accession(Primary)	P15056
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 84 kDa; Observed MW: 87 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P15056)

Gene Name BRAF

About Gene Symbol



Western blot analysis of B Raf in MDA-MB-468, T47D, MCF-7, HCT116, C2C12, 3T3 and HeLa lysates using B Raf antibody.

