

JAK1 (8B8) Mouse mAb


ABZ4242

Product Overview


Name	JAK1 (8B8) Mouse mAb
Accession(Primary)	P23458
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 133 kDa; Observed MW: 130 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20

Target information (P23458)


Gene Name	JAK1
About Gene Symbol	JAK1 (Janus Kinase 1) is a receptor tyrosine kinase involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of JAK1 (8B8) in C6, Ramos, Jurkat and HeLa lysates using Jak1 antibody.



Immunocytochemistry analysis of JAK1 (8B8) in HeLa using Jak1 antibody.



Immunoprecipitation analysis of JAK1 (8B8) in HeLa lysates using Jak1 antibody.

