

PKM2 (1A7) Mouse mAb


ABZ4243

Product Overview


Name	PKM2 (1A7) Mouse mAb
Accession(Primary)	P14618
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 58 kDa; Observed MW: 58 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (P14618)

Gene Name	PKM
About Gene Symbol	PKM (Pyruvate Kinase M1/2) is a receptor tyrosine kinase involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PKM2 in COS7, PC-12, C6, 3T3 and HeLa lysates using PKM2 antibody.



Immunocytochemistry analysis of PKM2 (1A7) in HeLa using PKM2 antibody.