

Placental Alkaline Phosphatase Mouse mAb


ABZ4245

Product Overview

Name	Placental Alkaline Phosphatase Mouse mAb
Accession(Primary)	P05187
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 58 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P05187)

Gene Name	ALPP
About Gene Symbol	ALPP (Alkaline Phosphatase, Placental) is a protein phosphatase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Placental Alkaline Phosphatase in HepG2, HeLa, CHO-K1 and CHO-K1 transfected pcDNA3.1HAPLAP lysates using Placental alkaline phosphatase (PLAP) antibody.