

DR5 Mouse mAb

ABZ4246


Product Overview

Name	DR5 Mouse mAb
Accession(Primary)	O14763
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Calculated MW	Calculated MW: 48 kDa; Observed MW: 40,48 kDa
Recommended Dilution	WB-1:500-1:1000


Target information (O14763)

Gene Name TNFRSF10B

About Gene Symbol



Western blot analysis of DR5 in SW480, K562, HCT116 and HeLa lysates using DR5 antibody.



Western blot analysis of TNFRSF10B (6F6) in DR5
in MCF-7 and T47D lysates using DR5 antibody.