


CD31 (5E2) Mouse mAb


ABZ4247

Product Overview

Name	CD31 (5E2) Mouse mAb
Accession(Primary)	P16284
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 83 kDa; Observed MW: 130 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (P16284)

Gene Name PECAM1**About Gene Symbol**



Immunohistochemistry analysis of paraffin-embedded Human tonsils using CD31/PECAM1 (5E2) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.