


HtrA2 / Omi Mouse mAb


ABZ4248

Product Overview

Name	HtrA2 / Omi Mouse mAb
Accession(Primary)	O43464
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 49 kDa; Observed MW: 37 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (O43464)

Gene Name HTRA2**About Gene Symbol**



Immunoprecipitation analysis of HtrA2 (8G11) in HeLa lysates using HtrA2/Omi antibody.