


NSE Mouse mAb

ABZ4251


Product Overview

Name	NSE Mouse mAb
Accession(Primary)	P09104
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:500-1:1000; ELISA-1:10000


Target information (P09104)

Gene Name ENO2**About Gene Symbol**

Western blot analysis of Enolase2 in SH-SY5Y and U87-MG lysates using Enolase2 antibody.



Western blot analysis of NSE (1H8) in rat Brain and mouse Brain lysates using Enolase2 antibody.



Observed Enolase2 levels in CHO-K1 cell lysates transfected with Enolase2 at different dilution.