

Cardiac Troponin I (1D7) Mouse mAb

ABZ4252


Product Overview

Name	Cardiac Troponin I (1D7) Mouse mAb
Accession(Primary)	P19429
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 24 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:500-1:1000; ELISA-1:10000


Target information (P19429)

Gene Name TNNI3

About Gene Symbol



sELISA standard Curve for Cardiac Troponin I:
 Capture antibody antibody to Cardiac Troponin I at
 2ug/ml and 2012951D7(1D7Dwas used for
 detecting.



Western blot analysis of Cardiac Troponin I (1D7)
in rat and mouse heart lysates using Cardiac
Troponin I (1D7) antibody