


Cardiac Troponin I (4D3) Mouse mAb

ABZ4254


Product Overview

Name	Cardiac Troponin I (4D3) Mouse mAb
Accession(Primary)	P19429
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 24 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:500-1:1000; ELISA-1:10000

Target information (P19429)

Gene Name TNNI3**About Gene Symbol**

Western blot analysis of Cardiac Troponin I (4D3) in rat and mouse heart lysates using Troponin I (2012954D3) antibody.



sELISA standard Curve for Cardiac Troponin I:
Capture antibody to Cardiac Troponin I at 2ug/ml
and 2012954D3(4D3F3) was used for detecting.