


Cardiac Troponin I (5C6) Mouse mAb

ABZ4255

Product Overview

Name	Cardiac Troponin I (5C6) Mouse mAb
Accession(Primary)	P19429
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 24 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:500-1:1000; ELISA-1:10000

Target information (P19429)

Gene Name TNNI3**About Gene Symbol**

Western blot analysis of Cardiac Troponin I (5C6) in rat and mouse heart lysates using Troponin I (2012955C6) antibody.

