


CaMKII alpha (5D9) Mouse mAb

ABZ4256

Product Overview

Name	CaMKII alpha (5D9) Mouse mAb
Accession(Primary)	Q9UQM7
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 54 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q9UQM7)

Gene Name CAMK2A**About Gene Symbol**

Western blot analysis of CaMKIIα in mouse Brain lysates using CaMKIIα antibody.

