

NSE Mouse mAb

ABZ4258


Product Overview

Name	NSE Mouse mAb
Accession(Primary)	P09104
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200


Target information (P09104)

Gene Name ENO2


About Gene Symbol




Western blot analysis of Enolase2 in COS7, C6 and HeLa lysates using Enolase2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human small cell carcinoma of lung tissue using NSE antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of NSE (4C7) in Enolase (4C7) in CHO-K1, rat brain, C6, 3T3, HeLa lysates using Enolase (4C7) antibody.



Immunocytochemistry analysis of NSE (4C7) (green) in HeLa using NSE (4C7) antibody, and DAPI (blue)

