

S100 alpha 6 (3E11) Mouse mAb

ABZ4262


Product Overview

Name	S100 alpha 6 (3E11) Mouse mAb
Accession(Primary)	P06703
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 10 kDa; Observed MW: 10 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P06703)

Gene Name S100A6

About Gene Symbol



Western blot analysis of S100A6 in COS7, MCF-7, HCT116 and A549 lysates using S100A6 antibody.

