

Placental Alkaline Phosphatase Mouse mAb


ABZ4263

Product Overview


Name	Placental Alkaline Phosphatase Mouse mAb
Accession(Primary)	P05187
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Calculated MW	Calculated MW: 58 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (P05187)


Gene Name	ALPP
About Gene Symbol	ALPP (Alkaline Phosphatase, Placental) is a protein phosphatase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of Placental Alkaline Phosphatase in HepG2 , A431 and MCF-7 lysates using Placental Alkaline Phosphatase antibody.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Placental Alkaline Phosphatase antibody, with DAB staining. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Placental Alkaline Phosphatase in HepG2 lysates using Placental Alkaline Phosphatase antibody



Western blot analysis of Placental Alkaline Phosphatase in HeLa, Human Embryo, HepG2 lysates using Placental Alkaline Phosphatase antibody.