

CD45 (8C1) Mouse mAb

ABZ4266


Product Overview

Name	CD45 (8C1) Mouse mAb
Accession(Primary)	P08575
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 145 kDa; Observed MW: 145 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200


Target information (P08575)

Gene Name Ptpnc


About Gene Symbol



Immunofluorescence analysis of CD45 (8C1) in Human livercancer tissue using CD45 (8C1) antibody(red),and DAPI (blue).



Immunohistochemical analysis of paraffin-embedded Human tonsils using CD45 (8C1) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Western blot analysis of CD45 (8C1) in Jurkat lysates using CD45 (8C1) antibody.