


MMP2 (3A8) Mouse mAb

ABZ4267


Product Overview

Name	MMP2 (3A8) Mouse mAb
Accession(Primary)	P08253
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 74 kDa; Observed MW: 74 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (P08253)

Gene Name MMP2**About Gene Symbol**

Western blot analysis of MMP2 (3A8) in Raw264.7 lysates using MMP2 antibody.



Immunofluorescence analysis of MMP2 (3A8) in Hela using MMP2 antibody, and DAPI (blue).