

Vimentin (4F8) Mouse mAb


ABZ4268

Product Overview


Name	Vimentin (4F8) Mouse mAb
Accession(Primary)	P08670
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human
Calculated MW	Calculated MW: 54 kDa; Observed MW: 57 kDa
Recommended Dilution	WB-1:5000-50000; IHC-1:500-1:4000; IF-1:100-1:500; IP-1:20

Target information (P08670)


Gene Name	VIM
About Gene Symbol	VIM (Vimentin) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Vimentin in Molt4, K562, COS7, Jurkat, HeLa and Vimentin negative (Daudi, Ramos, Raji) lysates using Vimentin antibody.



Immunocytochemistry analysis of Vimentin (4F8) in HeLa using Vimentin antibody.



Immunoprecipitation analysis of Vimentin (4F8) in HeLa lysates using Vimentin antibody.

