

E Cadherin (6B10) Mouse mAb

ABZ4269


Product Overview


Name	E Cadherin (6B10) Mouse mAb
Accession(Primary)	P12830
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 97 kDa; Observed MW: 135 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (P12830)


Gene Name CDH1

About Gene Symbol





Western blot analysis of E Cadherin (6B10) in mouse embryos and HCT116 lysates using ECadherin antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsils using ECadherin (6B1) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.