


PP1C alpha Mouse mAb

ABZ4272

Product Overview

Name	PP1C alpha Mouse mAb
Accession(Primary)	P62136
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200

Target information (P62136)

Gene Name PPP1CA**About Gene Symbol**

Western blot analysis of PPP1A in K562, 3T3 and HeLa lysates using PPP1A antibody.

