


## SMC1A Mouse mAb

ABZ4273

### Product Overview

|                             |                                              |
|-----------------------------|----------------------------------------------|
| <b>Name</b>                 | SMC1A Mouse mAb                              |
| <b>Accession(Primary)</b>   | Q14683                                       |
| <b>Host Species</b>         | Mouse                                        |
| <b>Clonality</b>            | Monoclonal                                   |
| <b>Isotype</b>              | IgG1                                         |
| <b>Species Reactivity</b>   | Human                                        |
| <b>Calculated MW</b>        | Calculated MW: 143 kDa; Observed MW: 143 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IF-1:50-1:200               |

### Target information (Q14683)

**Gene Name** SMC1A**About Gene Symbol**

Western blot analysis of SMC1A(Cterm) in K562, Jurkat, C6, 3T3 and HeLa lysates using SMC1A (Nterminus) antibody.

