


PP1C alpha Mouse mAb

ABZ4276

Product Overview

Name	PP1C alpha Mouse mAb
Accession(Primary)	P62136
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (P62136)

Gene Name PPP1CA**About Gene Symbol**

Western blot analysis of PPP1A in C6, COS7, K562, 3T3 and HeLa lysates using PPP1A antibody.

