


# Cytokeratin 18 Mouse mAb

ABZ4277

## Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | Cytokeratin 18 Mouse mAb                   |
| <b>Accession(Primary)</b>   | P05783                                     |
| <b>Host Species</b>         | Mouse                                      |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG2a                                      |
| <b>Species Reactivity</b>   | Human, Mouse, Cat                          |
| <b>Calculated MW</b>        | Calculated MW: 48 kDa; Observed MW: 48 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IF-1:50-1:200             |

## Target information (P05783)

**Gene Name** KRT18**About Gene Symbol**

Western blot analysis of Keratin 18 in COS7, 3T3 and HeLa lysates using Keratin 18 antibody.

