


CKMT1 Mouse mAb

ABZ4278

Product Overview

Name	CKMT1 Mouse mAb
Accession(Primary)	P12532
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P12532)

Gene Name CKMT1A**About Gene Symbol**

Western blot analysis of CKMT1 in rat Brain, 293T and CHO-K1 lysates using CKMT1 antibody.

