

Pyruvate Dehydrogenase E2 Mouse mAb


ABZ4279

Product Overview

Name	Pyruvate Dehydrogenase E2 Mouse mAb
Accession(Primary)	P10515
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Cat
Calculated MW	Calculated MW: 69 kDa; Observed MW: 69 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20

Target information (P10515)

Gene Name	DLAT
About Gene Symbol	DLAT (Dihydrolipoamide S-Acetyltransferase) is a metabolic dehydrogenase involved in adaptive immunity. It is an important research target in fibrotic disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Pyruvate Dehydrogenase E2 in COS7, Jurkat and HeLa lysates using Pyruvate Dehydrogenase E2 antibody.