

RAD9A Mouse mAb

ABZ4281


Product Overview

Name	RAD9A Mouse mAb
Accession(Primary)	Q99638
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (Q99638)

Gene Name RAD9A

About Gene Symbol



Western blot analysis of RAD9A in K562, 3T3 and HeLa lysates using RAD9A antibody.

