

# HIV 1 p24 (1F8) Mouse mAb

ABZ4294

## Product Overview

|                             |                           |
|-----------------------------|---------------------------|
| <b>Name</b>                 | HIV 1 p24 (1F8) Mouse mAb |
| <b>Accession(Primary)</b>   | P04591                    |
| <b>Host Species</b>         | Mouse                     |
| <b>Clonality</b>            | Monoclonal                |
| <b>Isotype</b>              | IgG1                      |
| <b>Species Reactivity</b>   | Virus                     |
| <b>Recommended Dilution</b> | ELISA-1:10000             |

## Target information (P04591)

|                          |     |
|--------------------------|-----|
| <b>Gene Name</b>         | gag |
| <b>About Gene Symbol</b> |     |