

NCAM1 Mouse mAb


ABZ4309

Product Overview


Name	NCAM1 Mouse mAb
Accession(Primary)	P13591
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Cat
Calculated MW	Calculated MW: 95kDa;Observed MW: 95kDa
Recommended Dilution	WB-1:1000-1:4000; IHC-1:100-1:400; FC-1:50-1:100; ELISA-1:10000

Target information (P13591)


Gene Name	NCAM1
About Gene Symbol	



Western blot analysis of NCAM1 in SH-SY5Y (1), COS-7 (2), and NIH3T3 (3) cell lysate using NCAM1 antibody.



Immunohistochemistry analysis of paraffin-embedded human breast cancer tissues using NCAM1 antibody.



Flow Cytometry analysis of Hela cells using NCAM1 antibody (green) and negative control (red).

