

Acetyl-Histone H3 (Lys9) Rabbit pAb


ABZ4317

Product Overview

Name	Acetyl-Histone H3 (Lys9) Rabbit pAb
Accession(Primary)	P68431
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 15 kDa; Observed MW: 15 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300; IF-1:50-1:200

Target information (P68431)

Gene Name	H3C1
About Gene Symbol	H3C1 (H3 Clustered Histone 1) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to el



Western blot analysis of Acetyl-Histone H3 (Lys9)
in COS7 and C6 lysates using Acetyl-Histone H3
(Lys9) antibody.