

Acetyl-Alpha Tubulin (Lys40) (6B5) Mouse mAb

ABZ4318


Product Overview

Name	Acetyl-Alpha Tubulin (Lys40) (6B5) Mouse mAb
Accession(Primary)	P68363
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100


Target information (P68363)

Gene Name TUBA1B


About Gene Symbol



Western blot analysis of Acetyl- α -Tubulin (Lys40) in rat spleen, Raw246.7 lysates using Acetyl- α -Tubulin (Lys40) (6B5) antibody.



Western blot analysis of Acetyl-Alpha Tubulin (Lys40) (6B5) in 293, C6 lysates using Acetyl- α -Tubulin (Lys40) (6B5) antibody.



Immunohistochemistry analysis of paraffin-embedded Human Brain Tissue using a tubulin(Acetyl Lys4) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.