

Acetyl-Histone H4 (Lys5) Rabbit mAb


ABZ4320

Product Overview


Name	Acetyl-Histone H4 (Lys5) Rabbit mAb
Accession(Primary)	P62805
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 11 kDa; Observed MW: 11 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:500-1:5000; IP-1:20-1:50

Target information (P62805)


Gene Name	H4C1
About Gene Symbol	




Immunocytochemistry analysis of Acetyl-Histone H4 (Lys5) (green) in HeLa using Acetyl-Histone H4 (Lys5) antibody, and DAPI (blue).



Western blot analysis of Acetyl-Histone H4 (Lys5) in 3T3 lysates using Acetyl-Histone H4 (Lys5) antibody. Observed band size: 11 kDa.



Western blot analysis of Acetyl-Histone H4 (Lys5) in A549, HL-60, U2OS, U87-MG lysates using Acetyl-Histone H4 (Lys5) antibody



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Acetyl-Histone H4 (Lys5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.