

Acetyl-p53 (Lys370) Rabbit mAb


ABZ4322

Product Overview


Name	Acetyl-p53 (Lys370) Rabbit mAb
Accession(Primary)	P04637
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 44 kDa; Observed MW: 53 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20

Target information (P04637)


Gene Name	TP53
About Gene Symbol	TP53 (Tumor Protein P53) is a tumor suppressor protein involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of Acetyl-p53 (Lys370) (green) in 3T3 using Acetyl-p53 (Lys370) antibody, and DAPI (blue).



Western blot analysis of Acetyl-p53 (Lys370) in HeLa, HeLa+TSA, 3T3, 3T3+TSA lysates using Acetyl-p53 (Lys370) antibody.



Western blot analysis of Acetyl-p53 (Lys370) in
rat Brain lysates using Acetyl-p53 (Lys370)
antibody