

Acetyl-p53 (Lys382) Rabbit pAb


ABZ4323

Product Overview

Name	Acetyl-p53 (Lys382) Rabbit pAb
Accession(Primary)	P04637
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 44 kDa; Observed MW: 53 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (P04637)

Gene Name	TP53
About Gene Symbol	TP53 (Tumor Protein P53) is a tumor suppressor protein involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Acetyl-p53 (Lys382) in
HepG2 lysates using Acetyl-p53 (Lys382) antibody.