

Acetyl-Histone H4 (Lys16) Rabbit pAb

ABZ4326


Product Overview

Name	Acetyl-Histone H4 (Lys16) Rabbit pAb
Accession(Primary)	P62805
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 11 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300; IF-1:50-1:200

Target information (P62805)

Gene Name H4C1

About Gene Symbol



Western blot analysis of Acetyl-Histone H4 (Lys16) in HeLa lysates untreated or treated with TSA lysates using Acetyl-Histone H4(Lys16) antibody(upper) and betaactin (5B7-) antibody(lower) TDY041.