

Acetyl-Histone H2B (Lys5) Rabbit pAb


ABZ4327

Product Overview


Name	Acetyl-Histone H2B (Lys5) Rabbit pAb
Accession(Primary)	P57053
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:10000

Target information (P57053)


Gene Name	H2BC3
About Gene Symbol	




Western blot analysis of Acetyl-Histone H2B (Lys5) in various cell lines lysates using Acetyl-Histone H2B-K5 antibody.



Immunofluorescence analysis of Acetyl-Histone H2B (Lys5) in C6 using Acetyl-Histone H2B-K5 antibody. C6 cells were treated by TSA, and DAPI (blue).



Immunofluorescence analysis of Acetyl-Histone H2B (Lys5) in HeLa cells using Acetyl-Histone H2B-K5 antibody. HeLa cells were treated by TSA, and DAPI (blue).



Immunofluorescence analysis of Acetyl-Histone H2B (Lys5) in NIH/3T3 cells using Acetyl-Histone H2B-K5 antibody, and DAPI (blue). NIH/3T3 cells were treated by TSA.