

Acetyl-Histone H2B (Lys15) Rabbit pAb


ABZ4328

Product Overview


Name	Acetyl-Histone H2B (Lys15) Rabbit pAb
Accession(Primary)	P57053
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:5000

Target information (P57053)


Gene Name	H2BC3
About Gene Symbol	




Western blot analysis of Acetyl-Histone H2B (Lys15) in various cell lines lysates using Acetyl-Histone H2BK15 antibody.




Immunohistochemistry analysis of paraffin-embedded rat lung using Acetyl-Histone H2B (Lys15) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.




Immunohistochemistry analysis of paraffin-embedded Human breast cancer using H2B K15ac antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human gastric cancer using H2B K15ac antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded mouse brain using H2B K15ac antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of Acetyl-Histone H2B (Lys15) in C6 using Acetyl-Histone H2BK15 antibody. C6 cells were treated by TSA, and DAPI (blue).