

Cleaved-Caspase 4/5 p20 Rabbit pAb

ABZ4332

Product Overview

Name	Cleaved-Caspase 4/5 p20 Rabbit pAb
Accession(Primary)	P49662
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 43 kDa; Observed MW: 47,22 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:40000

Target information (P49662)

Gene Name

CASP4

About Gene Symbol


