

Cleaved-Caspase 3 p12 Rabbit mAb


ABZ4333

Product Overview

Name	Cleaved-Caspase 3 p12 Rabbit mAb
Accession(Primary)	P42574
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 32 kDa; Observed MW: 32,12 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; ICC/IF-1:50-1:100; FC-1:100-1:200

Target information (P42574)

Gene Name	CASP3
About Gene Symbol	CASP3 (Caspase 3) is a cysteine protease (caspase) involved in apoptosis. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Caspase3 p12 in K562, 3T3, HeLa lysates using Cleaved-Caspase 3 p12 antibody.