

## Cleaved-Caspase 3 p17 Rabbit pAb


ABZ4335

### Product Overview


|                             |                                                |
|-----------------------------|------------------------------------------------|
| <b>Name</b>                 | Cleaved-Caspase 3 p17 Rabbit pAb               |
| <b>Accession(Primary)</b>   | P42574                                         |
| <b>Host Species</b>         | Rabbit                                         |
| <b>Clonality</b>            | Polyclonal                                     |
| <b>Isotype</b>              | IgG                                            |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                              |
| <b>Calculated MW</b>        | Calculated MW: 32 kDa; Observed MW: 17 kDa     |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; IHC-1:50-1:300; IF-1:50-1:200 |

### Target information (P42574)

|                          |                                                                                                                                                                                                                                                                   |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | CASP3                                                                                                                                                                                                                                                             |
| <b>About Gene Symbol</b> | CASP3 (Caspase 3) is a cysteine protease (caspase) involved in apoptosis. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of Cleaved-Caspase3 in various lysates using Cleaved-Caspase3 p17 antibody.



Immunohistochemistry analysis of paraffin-embedded mouse colon tissue using Cleaved-Caspase3 p17 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.