

Cleaved-Caspase 3 p17 Rabbit pAb


ABZ4336

Product Overview


Name	Cleaved-Caspase 3 p17 Rabbit pAb
Accession(Primary)	P42574
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 32 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P42574)

Gene Name	CASP3
About Gene Symbol	CASP3 (Caspase 3) is a cysteine protease (caspase) involved in apoptosis. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Cleaved-Caspase3 in Jurkat and 3T3 lysates using Cleaved-Caspase3 p17 antibody.



Western blot analysis of Cleaved-Caspase 3 p17 in Jurkat cells, untreated or treated with Staurosporine (1 μ M; 3 hr), using Cleaved-Caspase3 p17 antibody (upper) or β Actin (8F10) antibody(lower).