

Cleaved-Caspase 1 Rabbit pAb

ABZ4337


Product Overview

Name	Cleaved-Caspase 1 Rabbit pAb
Accession(Primary)	P29466
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 45 kDa; Observed MW: 25 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300; IF-1:50-1:200


Target information (P29466)

Gene Name CASP1

About Gene Symbol



Western blot analysis of Cleaved-Caspase 1 in NIH3T3 lysates using Cleaved-Caspase 1 antibody.



Immunohistochemical analysis of paraffin-embedded Human breast cancer using Cleaved-Caspase-1 antibody. High-pressure and temperature Sodium citrate pH 6.0 was used for antigen retrieval.