

Cleaved-PARP1 Rabbit mAb


ABZ4338

Product Overview

Name	Cleaved-PARP1 Rabbit mAb
Accession(Primary)	P09874
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 113 kDa; Observed MW: 89 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:200-1:500; FC-1:20-1:50; IP-1:20-1:50

Target information (P09874)

Gene Name	PARP1
About Gene Symbol	



Western blot analysis of Cleaved- PARP1 in K562
lysates using Cleaved-PARP1 antibody.