

Hydroxyl-Histone H2A (Tyr39) Rabbit mAb

ABZ4339


Product Overview


Name	Hydroxyl-Histone H2A (Tyr39) Rabbit mAb
Accession(Primary)	P04908
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:500-1:2000

Target information (P04908)


Gene Name H2AC4

About Gene Symbol





Western blot analysis of Hydroxyl-Histone H2A (Tyr39) in HeLa, A549, HL-60 lysates using Hydroxyl-Histone H2A (Tyr39) antibody. Observed band size: 14 kDa.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Hydroxyl-Histone H2A (Tyr39) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.