

Crotonyl-Histone H2B (Lys12) Rabbit mAb


ABZ4340

Product Overview

Name	Crotonyl-Histone H2B (Lys12) Rabbit mAb
Accession(Primary)	P33778
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:14 kDa;Observed MW: 14 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; ICC/IF-1:50-1:100

Target information (P33778)

Gene Name	H2BC4
About Gene Symbol	



Western blot analysis of Crotonyl-Histone H2B (Lys12) in HeLa, 3T3 lysates using Crotonyl-Histone H2B (Lys12) antibody.