

Phospho-CaMKII (Thr305) Rabbit pAb


ABZ4341

Product Overview


Name	Phospho-CaMKII (Thr305) Rabbit pAb
Accession(Primary)	Q9UQM7
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:5000

Target information (Q9UQM7)

Gene Name	CAMK2A
About Gene Symbol	



Western blot analysis of Phospho-CaMKII (Thr305) in various lysates using Phospho-CaMKII (Thr305) antibody.



Western blot analysis of Phospho-CaMKII (Thr305) in 453 lysates using Phospho-CaMKIIα/β/δ (T305) antibody.