

Phospho-HDAC6 (Ser22) Rabbit pAb


ABZ4344

Product Overview

Name	Phospho-HDAC6 (Ser22) Rabbit pAb
Accession(Primary)	Q9UBN7
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 131 kDa; Observed MW: 131 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:10000

Target information (Q9UBN7)

Gene Name	HDAC6
About Gene Symbol	HDAC6 (Histone Deacetylase 6) is a histone deacetylase involved in chromatin remodeling and epigenetics. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-HDAC6 (Ser22)
in 3T3 lysates using Phospho-HDAC6 (Ser22)
antibody.