

Phospho-PRKD2 (Ser876) Rabbit mAb

ABZ4346

Product Overview

Name	Phospho-PRKD2 (Ser876) Rabbit mAb
Accession(Primary)	Q9BZL6
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 97 kDa; Observed MW: 105 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (Q9BZL6)

Gene Name PRKD2

About Gene Symbol



