

Phospho-AP2M1 (Thr156) Rabbit mAb


ABZ4348

Product Overview

Name	Phospho-AP2M1 (Thr156) Rabbit mAb
Accession(Primary)	Q96CW1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q96CW1)

Gene Name	AP2M1
About Gene Symbol	AP2M1 (Adaptor Related Protein Complex 2 Subunit Mu 1) is a signaling adaptor protein involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cel



Western blot analysis of Phospho-AP2M1 (Thr156) in K562 lysates using Phospho-AP2M1 (Thr156) antibody.