

Phospho-FRS2 (Tyr436) Rabbit pAb

ABZ4349

Product Overview

Name	Phospho-FRS2 (Tyr436) Rabbit pAb
Accession(Primary)	Q8WU20
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Cat
Calculated MW	Calculated MW: 57 kDa; Observed MW: 65 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (Q8WU20)

Gene Name	FRS2
About Gene Symbol	

