

Phospho-TAOK1/2/3 (Ser181/Ser181/Ser177) Rabbit mAb


ABZ4350

Product Overview


Name	Phospho-TAOK1/2/3 (Ser181/Ser181/Ser177) Rabbit mAb
Accession(Primary)	Q7L7X3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 116 kDa; Observed MW: 105,116,138 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-1:50-1:100

Target information (Q7L7X3)


Gene Name	TAOK1/TAOK2/TAOK3
About Gene Symbol	



Western blot analysis of Phospho-TAOK1(Ser181)/TAOK2 (Ser181)/TAOK3 (Ser177) in 3T3 lysates using Phospho-TAOK1/2/3 (Ser181/Ser181/Ser177) antibody.



Western blot analysis of TAOK1 (Phospho-Ser181)/TAOK2 (Phospho-Ser181)/TAOK3 (Phospho-Ser177) in HeLa, A549, HL-60, U2OS, C6 lysates using TAOK1 (Phospho-Ser181)/TAOK2 (Phospho-Ser181)/TAOK3 (Phospho-Ser177) antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using TAOK1 (Phospho-Ser181)/TAOK2 (Phospho-Ser181)/TAOK3 (Phospho-Ser177) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.